

# Mapatumumab reference antibody



Catalog Number: 177143

## Product name

Mapatumumab reference antibody

## Specificity

Human TNFRSF10A / TRAILR1 / CD261

## Antibody description

Anti-TNFRSF10A / TRAILR1 / CD261 Antibody (Mapatumumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG2SA

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-TNFRSF10A / TRAILR1 / CD261 Antibody (mapatumumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-TNFRSF10A / TRAILR1 / CD261 Antibody (mapatumumab) is 98.89%, determined by SEC-HPLC.