

# Labetuzumab govitecan (Synonyms: ANTI-CEACAM5-SN 38)



Catalog Number: 177173

## Product name

Labetuzumab govitecan (Synonyms: ANTI-CEACAM5-SN 38)

## Specificity

Human CEACAM5 / CEA / CD66e

## Antibody description

Anti-CEACAM5 / CEA / CD66e Antibody (Labetuzumab govitecan)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




SEC-HPLC

The purity of Anti-CEACAM5 / CEA / CD66e Antibody (labetuzumab govitecan) is 95%, determined by SEC-HPLC.



HIC-HPLC

The drug-to-antibody ratio (DAR) of the Anti-CEACAM5 / CEA / CD66e Antibody (labetuzumab govitecan) is 8.0.



Bioactivity: FACS

Human CEACAM5 / CEA / CD66e HEK293 Cell Line were stained Anti-CEACAM5 / CEA / CD66e Antibody (labetuzumab govitecan) and negative control protein respectively, washed and then followed by PE and analyzed with FACS, The EC50 is 1.45 µg/mL.