

Guselkumab reference antibody



Catalog Number: 177178

Product name

Guselkumab reference antibody

Specificity

Human IL-23

Antibody description

Anti-IL-23 Antibody (Guselkumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-IL-23 Antibody (guselkumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-IL-23 Antibody (guselkumab) is 98.97%, determined by SEC-HPLC.




Bioactivity: ELISA

Immobilized human IL 23P40 His at 2 ug/mL can bind Anti-IL-23 Antibody (guselkumab), EC₅₀=0.009014 ug/mL.

Guselkumab reference antibody



Catalog Number: 177178



Function: Luciferase

Anti-IL-23 Antibody (guselkumab) induced report gene STAT3 activity was evaluated using Human IL23 HEK-Blue. The EC₅₀ was approximately 0.0095 ug/mL.