

Lifastuzumab reference antibody



Catalog Number: 177179

Product name

Lifastuzumab reference antibody

Specificity

Human NaPi2b / SLC34A2

Antibody description

Anti-NaPi2b / SLC34A2 Antibody (Lifastuzumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-NaPi2b / SLC34A2 Antibody (Lifastuzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-NaPi2b / SLC34A2 Antibody (Lifastuzumab) is 100%, determined by SEC-HPLC.