

Dezamizumab (Synonyms: ANTI - SAP MAB)



Catalog Number: 177186

Product name

Dezamizumab (Synonyms: ANTI - SAP MAB)

Specificity

Human Serum Amyloid P / SAP

Antibody description

Anti-Serum Amyloid P / SAP Antibody
(Dezamizumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

1×PBS, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG

Applications


ELISA

Validations



SDS-PAGE

Anti-Serum Amyloid P / SAP Antibody (dezamizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-Serum Amyloid P / SAP Antibody (dezamizumab) is more than 95%, determined by SEC-HPLC.