

# Zampilimab reference antibody



Catalog Number: 177188

## Product name

Zampilimab reference antibody

## Specificity

Human TGM2 / Transglutaminase 2

## Antibody description

Anti-TGM2 / Transglutaminase 2 Antibody  
(Zampilimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG4SP

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-TGM2 / Transglutaminase 2 Antibody (zampilimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-TGM2 / Transglutaminase 2 Antibody (zampilimab) is more than 95%, determined by SEC-HPLC.