

# Xentuzumab reference antibody



Catalog Number: 177193

## Product name

Xentuzumab reference antibody

## Specificity

Human IGF-1

## Antibody description

Anti-IGF-1 Antibody (Xentuzumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-IGF-1 Antibody (xentuzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95.3%.



## SEC-HPLC

The purity of Anti-IGF-1 Antibody (xentuzumab) is 98.72%, determined by SEC-HPLC.



## Bioactivity: ELISA

Immobilized human IGF I Protein, His Tag at 2 ug/mL can bind Anti-IGF-1 Antibody (xentuzumab), EC50=0.009806 ug/mL.