

# Donanemab reference antibody



Catalog Number: 177196

## Product name

Donanemab reference antibody

## Specificity

Human Amyloid Beta

## Antibody description

Anti-Amyloid Beta Antibody (Donanemab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications

ELISA

## Validations



## SDS-PAGE

Anti-Amyloid Beta Antibody (donanemab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.