

# Donanemab reference antibody



Catalog Number: 177196

## Product name

Monoclonal

Donanemab reference antibody

## Ig Type

## Specificity

IgG1

Human Amyloid Beta

## Applications

## Antibody description

ELISA

Anti-Amyloid Beta Antibody (Donanemab)

## Validations

## Preparation

Recombinant expression and purified from CHO cells.


## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality



## SDS-PAGE

Anti-Amyloid Beta Antibody (donanemab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.