

# Ecromeximab reference antibody



Catalog Number: 177210

## Product name

Ecromeximab reference antibody

## Specificity

Human Disialoganglioside GD3

## Antibody description

Anti-Disialoganglioside GD3 Antibody  
(Ecromeximab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-GD3 Antibody (ecromeximab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-GD3 Antibody (ecromeximab) is 96.37%, determined by SEC-HPLC.