

Etaracizumab reference antibody



Catalog Number: 177221

Product name

Etaracizumab reference antibody

Specificity

Human Integrin α V β 3 (ITGAV & ITGB3)

Antibody description

Anti-Integrin α V β 3 (ITGAV & ITGB3) Antibody (Etaracizumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

1×PBS, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG

Applications


ELISA

Validations



SDS-PAGE

Anti-Integrin α V β 3 (ITGAV & ITGB3) Antibody (etaracizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-Integrin α V β 3 (ITGAV & ITGB3) Antibody (etaracizumab) is more than 95%, determined by SEC-HPLC.