

# Pritoxaximab reference antibody



Catalog Number: 177233

## Product name

Pritoxaximab reference antibody

## Specificity

E.coli Shiga toxin (E.coli)

## Antibody description

Anti-Shiga toxin (E.coli) Antibody (Pritoxaximab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



### SEC-HPLC

The purity of Anti-Shiga toxin (E.coli) Antibody (pritoxaximab) is 100.00%, determined by SEC-



### HPLC.

### SDS-PAGE

Anti-Shiga toxin (E.coli) Antibody (pritoxaximab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.