

Pritoxaximab reference antibody



Catalog Number: 177233

Product name

Pritoxaximab reference antibody

Specificity

E.coli Shiga toxin (E.coli)

Antibody description

Anti-Shiga toxin (E.coli) Antibody (Pritoxaximab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SEC-HPLC

The purity of Anti-Shiga toxin (E.coli) Antibody (pritoxaximab) is 100.00%, determined by SEC-



HPLC.

SDS-PAGE

Anti-Shiga toxin (E.coli) Antibody (pritoxaximab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.