

Anti-LTA antibody

Catalog Number: 175007

Product name

Polyclonal

Anti-LTA antibody

Ig Type**Specificity**

Rabbit IgG

Human, Rat

Applications**Antibody description**

WB

Rabbit polyclonal antibody to LTA

Dilutions**Preparation**

WB: 1:1,000

This antigen of this antibody was synthetic peptide within human lt-alpha aa 51-104.

Validations**Formulation**

Liquid, 1*PBS (pH7.4), 0.2% BSA, 25% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Fig1: Western blot analysis on LT-alpha (TNF- β) using anti- LT-alpha polyclonal antibody