

ZNF586 antibody



Catalog Number: 117007

Product name

ZNF586 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF586 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, COLO 320 cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of ZNF586 recombinant protein (Accession Number: NM_001077426). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117007(ZNF586 Antibody) at dilution of 1:600