

Mouse Tgm2 (His tag) recombinant protein



Catalog Number: 518007

General Information

Gene Name Synonym

Protein-glutamine gamma-glutamyltransferase 2, Tgm2, Tissue transglutaminase, Transglutaminase C, TGase-2

Protein Construction

A DNA sequence encoding the Mouse Tgm2 Ala2-Ala686 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 20% Glycerol, pH7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Tgm2 (His tag) recombinant protein