

CHD1 Polyclonal Antibody



Catalog Number: 162407

Product name

CHD1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CHD1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1501-1710 of human CHD1 (NP_001261.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - CHD1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using CHD1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.