

nucleobindin 1 antibody



Catalog Number: 113407

Product name

nucleobindin 1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:50-1:500

Validations

Antibody description

nucleobindin 1 Rabbit Polyclonal antibody.
Positive WB detected in HepG2 cells, MCF7 cells.
Positive IF detected in HepG2 cells. Positive IHC detected in human spleen tissue, human kidney tissue, human liver tissue, human placenta tissue.
Observed molecular weight by Western-blot: 63 kDa



Preparation

This antibody was obtained by immunization of nucleobindin 1 recombinant protein (Accession Number: BC015787). Purification method: Antigen affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113407(NUCB1 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:




Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:113407(NUCB1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

nucleobindin 1 antibody

Catalog Number: 113407



Immunohistochemical of paraffin-embedded human spleen using Catalog No:113407(NUCB1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human spleen using Catalog No:113407(NUCB1 antibody) at dilution of 1:100 (under 40x lens)