

Anti-PAOX antibody

Catalog Number: 176407

Product name

Anti-PAOX antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PAOX

Preparation

This antigen of this antibody was synthesized peptide derived from paox . at aa range: 260-340

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-300

Validations




Fig1:: Western Blot analysis of extracts from rat stomach, K562 cells, using PAOX Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100