

Anti-ZBTB45 antibody

Catalog Number: 175807

Product name

Anti-ZBTB45 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to ZBTB45

Preparation

This antigen of this antibody was recombinant protein of zbtb45

Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:1000

IHC-P: 1:200-1:500

Validations




Fig1.; Western blot analysis of 1) HeLa, 2) Jurkat, 3) Mouse Brain, 4) Mouse Kidney, 5) Rat Brain using ZBTB45 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at




Fig2.; Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using ZBTB45 Polyclonal Antibody.