

# Anti-SNX9 Rabbit antibody



Catalog Number: 166807

## Product name

Anti-SNX9 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

SNX9 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human SH3PX1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/20

FC: 1/20

IP: 1/20

## Validations



Western blot detection of SH3PX1 in K562, HeLa cell lysates using SH3PX1 Rabbit pAb(1:1000 diluted). Predicted band size: 67kDa. Observed band size: 67kDa.