

Anti-Prealbumin Rabbit antibody



Catalog Number: 168807

Product name

Anti-Prealbumin Rabbit antibody

Specificity

Human

Antibody description

Prealbumin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Prealbumin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/1000

FC: 1/20

Validations



Western blot detection of Prealbumin in Rat Brain lysates using Prealbumin Rabbit mAb(1:1000 diluted).Predicted band size:16kDa.Observed band size:16kDa.