

Human FUT8 (aa 68-575, His Tag) recombinant protein



Catalog Number: 502107

General Information

Gene Name Synonym

Fucosyltransferase 8; GDP-L-Fuc:N-acetyl-beta-D-glucosaminide alpha1,6-fucosyltransferase; GDP-fucose--glycoprotein fucosyltransferase; Glycoprotein 6-alpha-L-fucosyltransferase

Protein Construction

A DNA sequence encoding the human FUT8 isoform 1 (Q9BYC5-1) (Arg 68-Lys 575) was fused with a polyhistidine tag at the carboxy-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Measured by its ability to hydrolyze the donor substrate GDP fucose.
The specific activity is >0.75 pmoles/min/ μg .

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Arg 68

Molecular Mass

The recombinant human FUT8 consists of 518 amino acids and has a calculated molecular mass of 60 kDa. It migrates as an approximately 55 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 8.0, 10% gly

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human FUT8 Protein (aa 68-575, His Tag) SDS-PAGE