

# KLHL4 antibody

Catalog Number: 112107

## Product name

KLHL4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

KLHL4 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 75-80kd

## Preparation

This antibody was obtained by immunization of KLHL4 recombinant protein (Accession Number: NM\_019117). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112107(KLHL4 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human brain using Catalog No:112107(KLHL4 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:112107(KLHL4 antibody) at dilution of 1:50 (under 40x lens)