

# NEK9 antibody

Catalog Number: 113107

---

## Product name

NEK9 antibody

## Specificity

Human; other species not tested.

## Antibody description

NEK9 Rabbit Polyclonal antibody. Positive IF detected in COS-7 cells.

## Preparation

This antibody was obtained by immunization of NEK9 recombinant protein (Accession Number: BC009336). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:20-1:200

## Validations



Immunofluorescent analysis of (10% Formaldehyde) fixed COS-7 cells using Catalog No:113107(NEK9 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)