

Anti-Atp1a1 antibody



Catalog Number: 176107

Product name

Anti-Atp1a1 antibody

Specificity

Rat

Antibody description

Rabbit polyclonal antibody to Atp1a1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from rat atp1 alpha1/na+k+ atpase1. aa range:15-64

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western Blot analysis of HeLa cells using Na⁺/K⁺-ATPase α1 Polyclonal Antibody