

Anti-NFκB p105/p50 (5E3) Mouse antibody



Catalog Number: 168107

Product name

Anti-NFκB p105/p50 (5E3) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

NFκB p105/p50 (5E3) Mouse monoclonal antibody

Dilutions

WB: 1/500

Preparation

Antigen: Recombinant human NF-κB1 p105/p50 protein.

Validations

Formulation


PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of NF-κB1 p105/p50 in K562, C6, 3T3 and HeLa cell lysates using NF-κB1 p105/p50 mouse mAb(dilution 1:500). Predicted band size: 120, 50kDa. Observed band size: 120, 50kDa.