

# Phospho-ATR-S428 pAb



Catalog Number: 166507

---

**Product name**

Rabbit IgG

Phospho-ATR-S428 pAb

**Applications****Specificity**

WB, IF

Human

**Dilutions****Antibody description**

WB 1:500 - 1:2000

Polyclonal antibody to Phospho-ATR-S428


IF 1:50 - 1:200

**Preparation**

Antigen: A synthetic phosphorylated peptide around S428 of human ATR (NP\_001175.2).

**Validations****Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Immunofluorescence - Phospho-ATR-S428 pAb

**Clonality**

Polyclonal

Immunofluorescence analysis of U2OS cells using Phospho-ATR-S428 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

**Ig Type**