

RND3 antibody



Catalog Number: 107507

Product name

RND3 antibody

Specificity

Human, Pig; other species not tested.

Antibody description

RND3 Mouse Monoclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of RND3 recombinant protein (Accession Number: NM_005168). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107507(RND3 Antibody) at dilution of 1:2000