

# Anti-KAT13D/CLOCK Mouse antibody



Catalog Number: 168507

## Product name

Anti-KAT13D/CLOCK Mouse antibody

## Specificity

Human

## Antibody description

KAT13D/CLOCK Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human CLOCK expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ICC/IF

## Dilutions


ICC: 1/200 - 1/1000

ELISA: 1/10000

## Validations



Immunofluorescence analysis of HeLa cells using CLOCK mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Western blot analysis using CLOCK mouse mAb against CLOCK(AA: 200-465)-hIgGfc transfected HEK293 cell lysate.