

Human HPCAL1/BDR1 (His tag) recombinant protein



Catalog Number: 518507

General Information

Gene Name Synonym

Hippocalcin-Like Protein 1, Calcium-Binding Protein BDR-1, HLP2, Visinin-Like Protein 3, VILIP-3, HPCAL1, BDR1

Protein Construction

A DNA sequence encoding the Human HPCAL1/BDR1 Met1-Phe193 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human HPCAL1/BDR1 (His tag) recombinant protein