

# EFHA1 antibody



Catalog Number: 110207

## Product name

EFHA1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

EFHA1 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in mouse pancreas tissue. Observed molecular weight by Western-blot: 50kd

## Preparation

This antibody was obtained by immunization of EFHA1 recombinant protein (Accession Number: NM\_152726). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:110207(EFHA1 antibody) at dilution of 1:800



Immunofluorescent analysis of HepG2 cells, using EFHA1 antibody Catalog No:110207 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).