

LGI2 antibody



Catalog Number: 112207

Product name

LGI2 antibody

Immunogen

[Human LGI2 Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

LGI2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 62kd

Preparation

This antibody was obtained by immunization of LGI2 recombinant protein (Accession Number: NM_018176). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112207(LGI2 antibody) at dilution of 1:800