

# Anti-FLT-3/CD135/FLK-2 antibody



Catalog Number: 106207

## Product name

Anti-FLT-3/CD135/FLK-2 antibody

## Immunogen

[Human FLT-3/CD135/FLK-2 \(His Tag\) recombinant protein](#)

## Specificity

Human FLT3 / CD135 / FLK2

**No cross-reactivity** in ELISA with E.coli cell lysate

## Antibody description

Mouse monoclonal to FLT-3/CD135/FLK-2

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human FLT3 / CD135 / FLK2 (rh FLT3 / CD135 / FLK2; NP\_004110.2; Met 1-Asn 541) and conjugated with APC under optimum conditions, the unreacted APC was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of

receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

FCM

## Dilutions

FCM

## Validations



CRABP2 Antibody, Mouse MAb,  
Immunohistochemistry

FLT3 / CD135 / FLK2 Antibody (APC), Mouse  
MAb, Flow cytometric analysis