

# FGF12 antibody



Catalog Number: 107207

## Product name

FGF12 antibody

## Specificity

Human; other species not tested.

## Antibody description

FGF12 Mouse Monoclonal antibody. Positive WB detected in Human brain, human brain tissue. Observed molecular weight by Western-blot: 20 kDa

## Preparation

This antibody was obtained by immunization of FGF12 recombinant protein (Accession Number: NM\_004113). Purification method: Caprylic Acid/Ammonium Sulfate Precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107207(FGF12 antibody) at dilution of 1:1000