

Anti-Placental alkaline phosphatase (PLAP) Mouse antibody



Catalog Number: 168207

Product name

Anti-Placental alkaline phosphatase (PLAP) Mouse antibody

Specificity

Human

Antibody description

Placental alkaline phosphatase (PLAP) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant Placental alkaline phosphatase (PLAP) protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot analysis of extracts from HepG2, Hela, CHO-K1 and CHO-K1 transfected pcDNA3.1-HA-PLAP cell lysates using Placental alkaline phosphatase (PLAP) mouse mAb (1:1000 diluted). Predicted band size: 70KDa. Observed band size: 70KDa.