

Anti-RBP4 antibody



Catalog Number: 104607

Product name

Monoclonal

Anti-RBP4 antibody

Ig Type

Rabbit IgG

Immunogen

[Mouse RBP4 \(His Tag\) recombinant protein](#)

Applications

FCM

Specificity

Dilutions

Mouse RBP4

FCM

Antibody description

Validations

Rabbit monoclonal to RBP4

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse RBP4 (rM RBP4; NP_035385.1; Met1-Leu201) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.


Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality



RBP4 Antibody (FITC), Rabbit MAb, Flow cytometric

Flow cytometric analysis of Mouse RBP4 expression in spleen lymphocytes. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), and then stained with FITC Rabbit anti-Mouse RBP4. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.