

# EIF2AK3 Monoclonal Antibody



Catalog Number: 166607

## Product name

EIF2AK3 Monoclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Monoclonal antibody to EIF2AK3

## Preparation

Antigen: Recombinant protein of human EIF2AK3

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C (regular) or -20°C (long term). Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - EIF2AK3 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using EIF2AK3 antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.