

# c-SRC antibody



Catalog Number: 109607

## Product name

c-SRC antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

c-SRC Rabbit Polyclonal antibody. Positive WB detected in A431 cells, CHO cells, COLO 320 cells, Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of c-SRC recombinant protein (Accession Number: XM\_047440405). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109607(c-SRC Antibody) at dilution of 1:1000