

SLC18A1 antibody



Catalog Number: 115307

Product name

SLC18A1 antibody

Immunogen

[Human SLC18A1 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

SLC18A1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A549 cells, C6 cells, HeLa cells, HepG2 cells. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 50 kDa, 56 kDa

Preparation

This antibody was obtained by immunization of SLC18A1 recombinant protein (Accession Number: BC009387). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115307(SLC18A1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human placenta using Catalog No:115307(SLC18A1 antibody) at dilution of 1:100 (under 10x lens)



SLC18A1 antibody



Catalog Number: 115307

Immunohistochemical of paraffin-embedded human placenta using Catalog

No:115307(SLC18A1 antibody) at dilution of 1:100 (under 40x lens)