

TRIM25 antibody



Catalog Number: 116307

Product name

TRIM25 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRIM25 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in MCF7 cells, Jurkat cells, rat spleen tissue. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 71 kDa

Preparation

This antibody was obtained by immunization of TRIM25 recombinant protein (Accession Number: NM_005082). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116307(TRIM25 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using TRIM25 antibody Catalog No:116307 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



DNA dye).

TRIM25 antibody



Catalog Number: 116307

IP Result of anti-TRIM25 (IP:Catalog No:116307,

4ug; Detection:Catalog No:116307 1:600) with
K-562 cells lysate 3200ug.