

# Anti-LDHA Rabbit antibody



Catalog Number: 167307

## Product name

Anti-LDHA Rabbit antibody

WB: 1/1000

IHC: 1/50

## Specificity

Human, Mouse, Rat

ICC/IF: 1/50

FC: 1/20

IP: 1/20

## Antibody description

LDHA Rabbit polyclonal antibody

## Validations

## Preparation

Antigen: A synthetic peptide of human Lactate Dehydrogenase

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal


## Ig Type

Rabbit IgG


## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions



Western blot detection of LDHA in K562, Rat Brain, C6, 3T3, HeLa cell lysates using LDHA Rabbit pAb(1:1000 diluted). Predicted band size: 37kDa. Observed band size: 37kDa.



Immunocytochemistry of Lactate Dehydrogenase(green) in HeLa cells using Lactate Dehydrogenase Rabbit pAb at dilution 1/50, and DAPI(blue)