

Anti-Neurologin 1/NLGN1 antibody



Catalog Number: 101707

Product name

Anti-Neurologin 1/NLGN1 antibody

Immunogen

[Human Neurologin 1/NLGN1 \(His Tag\) recombinant protein](#)

Specificity

Human Neurologin 1 / NLGN1

Antibody description

Rabbit monoclonal to Neurologin 1/NLGN1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Neurologin 1 / NLGN1 (rh Neurologin 1 / NLGN1; NP_055747.1; Met1-Ser677).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 1-10 µg/mL

Validations



Neurologin 1 / NLGN1 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human NLGN1 in human brain with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).