

# RIMS4-Specific antibody



Catalog Number: 114707

## Product name

RIMS4-Specific antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RIMS4-Specific Rabbit Polyclonal antibody.  
Positive WB detected in mouse brain tissue.  
Observed molecular weight by Western-blot: 29kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_182970).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114707(RIMS4-Specific antibody) at dilution of 1:1000