

# ERCC4 Polyclonal Antibody



Catalog Number: 165707

## Product name

ERCC4 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ERCC4

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 600-700 of human ERCC4 (NP\_005227.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB 1:500 - 1:2000


IP 1:20 - 1:50

## Validations



## Western blot - ERCC4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ERCC4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 30s.

## Immunoprecipitation - ERCC4 Polyclonal Antibody

Immunoprecipitation analysis of 150ug extracts of 293T cells using 3ug ERCC4 antibody . Western blot was performed from the immunoprecipitate using ERCC4 antibody at a dilution of 1:1000.