

# Anti-Nectin 2 Rabbit antibody



Catalog Number: 167707

## Product name

Anti-Nectin 2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Nectin 2 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Nectin 2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/50

FC: 1/20

IP: 1/20-1/100

## Validations



Western blot detection of Nectin 2 in HeLa, A549, HL-60, U251, U87-MG cell lysates using Nectin 2 Rabbit mAb(1:1000 diluted). Predicted band size: 58kDa. Observed band size: 70-80kDa.