

Human NECAP2 (His tag) recombinant protein



Catalog Number: 519707

General Information

Gene Name Synonym

Adaptin ear-binding coat-associated protein 2,
NECAP endocytosis associated protein 2, NECAP2

Protein Construction

A DNA sequence encoding the Human NECAP2
Met1-Phe263 is expressed with a 6His tag at the
N-terminus, 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM PB \square
150mM NaCl,pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human NECAP2 (His tag) recombinant protein