

Anti-TSNAX antibody



Catalog Number: 143017

Product name

Rabbit IgG

Anti-TSNAX antibody

Applications

Specificity

IHC-P, ICC/IF, IF

Human TSNAX

Dilutions

Antibody description

IHC-P: 0.1-2µg/mL

Rabbit Polyclonal to Human TSNAX

ICC/IF: 1-5µg/mL

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human TSNAX, and purified by antigen affinity chromatography.

Validations

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.


Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.


Clonality

Polyclonal

Ig Type



Human TSNAX Immunohistochemistry(IHC)



15658

Human TSNAX Immunofluorescence(IF) 15659