

# NAGS antibody

Catalog Number: 113017

## Product name

NAGS antibody

IHC: 1:20-1:200

IF: 1:10-1:100

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NAGS Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive IHC detected in human liver. Positive WB detected in HepG2 cells, Jurkat cells, L02 cells. Observed molecular weight by Western-blot: 58-66kd

## Preparation

This antibody was obtained by immunization of NAGS recombinant protein (Accession Number: XM\_011524439). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113017(NAGS antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:113017(NAGS antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human liver using Catalog No:113017(NAGS antibody) at dilution of 1:50 (under 10x lens)

# NAGS antibody

Catalog Number: 113017

---

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using NAGS antibody Catalog No:113017 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).