

TRIM63 Polyclonal Antibody



Catalog Number: 164017

Product name

TRIM63 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TRIM63

Preparation

Antigen: A synthetic peptide of human TRIM63

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - TRIM63 Polyclonal Antibody

Western blot analysis of extracts of mouse heart, using TRIM63 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunofluorescence - TRIM63 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using TRIM63 antibody. Blue: DAPI for nuclear staining.